

Reanibex 100 with connectivity Bexen Cardio

The Reanibex 100 is a portable, user-friendly automated external defibrillator (AED) which empowers non-professionals to provide fast, effective first aid treatment to cardiac arrest victims, adults or children.

Choose between WiFi or Sigfox connectivity, with a free Sigfox license included for the first 5 years.



Advantages of the Reanibex 100 with connectivity

Reanibex Data Net

Synchronize your AED with the Reanibex Data Net. Here, you can customize the settings of your AED, including the volume, CPR sequence, energy levels, and patient type. Plus, you can effortlessly check the battery status and the expiration date of the electrodes. Use of the data network is free of charge

/ Metronome

Feedback on the pace of CPR. The AED helps you if you do not resuscitate at the rhythm of the metronome. chances are high that due to this feature, the rescuer will resuscitate at a better pace. There is a high chance of a better outcome for the victim.

★ Child button

This button turns the AED into an AED for children under the age of 8. With the same electrodes a child can be helped without special child electrodes.

Water- and dust-resistant

The high IP rating of 55, guarantees resistance to water and dust. This makes this AED suitable in humid and / or dusty environments.

Semi or fully-automatic

In the semi-automatic version the rescuer has to push the shock button to deliver the shock. In the Fully-Automatic version the shock is delivered automatically when needed.

Durability

Engineered with robust materials and built to withstand challenging environments, the Reanibex 100 AED is designed for long-term reliability.

When do you receive alerts?

- Battery issues: When the battery is low or the battery needs to be replaced.
- Electrode Issues: When the electrodes are used or (almost) expire.
- The device cover is left open: In this case, the power consumption is higher.
- The device has been used with a patient: (or a patient is connected).
- **Shockable rhythm**: When the device detects a shockable rhythm.
- **Errors** that affect the functionality of the device.
- When there is a warning regarding the operation of the device.





General

Dimensions Without handle 217 mm (Width) x 235 mm (Depth) x 96 mm (High)

With handle 217 mm (Width) x 285 mm (Depth) x 96 mm (High)

Weight Equipment without battery: 1,7 kg Equipment without handle, with

disposable pads and battery: < 2kg Equipment with handle, with

disposable pads and battery: < 2,15 kg

AED Defibrillator

Operating mode Semi or fully automatic

Waveform Output energy Biphasic truncated exponential adapted to patient's impedance

Charging Time Pads 150 - 200 Joules for adult patients and 50 - 90 Joules for pediatric patients.

types Defibrillable 0 seconds - Precharged

rhythms CPR Help Disposable pre-connected pads for adult and pediatric. Resuscitation guidelines Ventricular fibrillations and rapid ventricular tachycardias. Metronome with compressions rate feedback in real time.

Factory set Guidelines 2015 (ERC/AHA) and its review of 2017.

Data storage

Internal memory 10 episodes of at least 15 hours of ECG signal and events.

Data review Web application Reanibex Data Net for downloading, reproducing, storing and

reviewing recorded episodes (optional).

Applications

USB Sync PC application to synchronize the device with the Reanibex Data Net application. NFC Reader

NFC application for Android devices to check the information and status of

the device and its accessories.

Web application to download, display, manage, store and analyze episodes. Reanibex Data Net

Communications

Wi-Fi 802.11 b/g/n Security Protocols: WPA2/WPA Personal, WPA2/WPA Enterprise

and WEP (for alerts and geolocation).

Frequency: between 868 and 868.2 MHz, Transmission power (typ): Between Sigfox

12.5 and 13.5 dBm (for alerts and geolocation).

USB Port USB 2.0 (for equipment upgrade and episodes sending).

Battery

Type Disposable battery LiMnO2 4.2 A/h, 12 V

Capacity More than 350 shocks at 200 J or more than 36 hours of monitoring at 25 °C.

Shelf life 5 years

Capacity At least 60 shocks 200J

Placement of electrodes Anterior/lateral, for children anterior/posterior

Electrode surface 100 cm2 Cable length 1 meter

Environment

Operating Temperature From 0 to 50 °C Storage temperature From -25 to 70 °C

Humidity 0 to 95% Altitude 0 to 5.000 m Resistance to solids/liquids IP55

EMC EN 60601-1-11:2015 and EN 60601-1-2:2015

Vibrations / Shocks EN 60601-1-11:2015

